
## RespiVet NANO II Respiratory Monitor

**Can be connected to a ventilator or directly to a Endotracheal Tube.**

### **Simple operation**

As long as the Respiratory Monitor is connected with the endotracheal tube and anesthesia machine, everything is ready.

## **Responsive**

Once the pet stops breathing for more than 15 seconds, the Respiratory Monitor sends an alarm. If the pet is breathing again during the 15 seconds, it can automatically wake up the device to continue to monitor.

## **Auto Start**

After connecting with the pet, please activate the upper button to start the device. After starting, the red light and green light flicker once, and an alarm sounds; During the monitoring process, the green light flashes once per breath, and an alarm sounds.

## **Auto Close**

After disconnection from the pet or if the pet stops breathing, the Respiratory Monitor will pause for 15 seconds, and then automatically close after the alertor sounds for 20 seconds. During the alarm, the red light flashes continuously.

## **Built in lithium battery**

The Respiratory Monitor uses lithium battery, the battery cannot be replaced or repaired.

## **Easy to clean**

After use, wipe with alcohol, and then dry with cotton cloth. Do not soak or disinfect with an pressure sterilizer.

### **Storage in closed areas**

When not in use, please keep it in a plastic bag or drawer to prevent the flowing gas from triggering the sensitive alarm function.

### **Convenient charging**

When the device is connected to USB for charging, first enter the start mode. The red and green lights flicker and the alarm sounds, and then enter the charging mode. Red light when charging and green light when fully charged.

### **Warning:**

The device will not automatically wake up if breathing again when the equipment alarms! If continue to use, wait for the alarm to end and start again to use it! The equipment during the use or alarm process cannot be manually switched off! The breathing condition cannot be monitored when the device is turned off, unless the start key is touched!

### **Configuration Attachment:**

USB (5 V): 1 pcs


User Manual: 1 pcs

**One year warranty period.**

## **Why use the Respiratory Monitor?**

Almost all Pet anesthetics inhibit respiration in a dose related manner, which makes the Respiratory Monitor an ideal tool for detecting respiration and evaluating the depth of anesthesia. Due to the lack of detection, most emergencies during anesthesia can be attributed to human error. The Respiratory Monitor is always on alert, which can not only prevent accidental death, but also better manage the use of anesthetic.

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