

Item No.	Autoclavable	Temperature [°C]	max. Temperature [°C]	Duration [min]	max. Duration [min]	min. Drying Time [min]
305004	✗					
306200	✗					
306201	✗					
30495000	✓	134		4	30	
30495031	✓	134		4	30	
30500005	✓	134		4	30	
30495033	✓	134		4	30	
305002	✓	134		4	30	
323135	✗					
323106	✓	134		3		15
323139	✗					
30495013	✓	132	137	3	18	10
30495014	✗					
30495050	✓	132	138	4	20	15
30495016	✓	132	137	5		20
305300	✗					
30526001	✗					