

# HAEMATOKRIT 200 TABLETOP CENTRIFUGE

## SEVEN MINUTES OF SEPARATION

**Hettich**  
CENTRIFUGES

### Compact tabletop centrifuge

Spin 24 capillaries per run in a compact frame with all-metal containment.

### Evaluation Lid

The lid doubles as an evaluation disc allowing you to determine hematocrit values without having to remove the capillaries disturbing the sedimentation.

### Run-Silent™ Operation

Excessive noise in the laboratory can lead to distractions, stress, anxiety, high blood pressure and chronic fatigue. Our centrifuges are as quiet as a normal conversation and will not increase the noise levels within your lab.

### Easy-Lift™

Our Easy-Lift Lid, featuring stay-in-place positioning and lid drop protection, is ergonomically designed to limit repetitive stress injuries, spillage and accidental sample disruption when loading and unloading your centrifuge.

### Safety-Lock™

Our Safety-Lock lids are self-closing, auto-locking and come with either single or dual steel lid locks. This ensures a safe, tight closure of the centrifuge chamber during operation. The centrifuge will not operate if lid is not secure. For added safety, our centrifuges feature **all-metal containment** and adhere to all national and international safety standards and certifications.

### Five (5) Year Warranty

FDA registered centrifuge manufacturer



MAX CAPACITY:  
24 x Capillaries



MAX RPM/RCF:  
13,000/16,060



DIMS. (HxWxD) :  
9.0" x 10.3" 13.9"


→ [www.hettweb.com/haematokrit200](http://www.hettweb.com/haematokrit200)

# HETTICH HAEMATOKRIT 200


## APPLICATION PACKAGES

All Hettich centrifuge packages come complete with everything you need to perform your application – centrifuge, rotor, buckets, adapters, inserts, bio-containment lids (optional where applicable), five (5) year warranty and free support.

### HAEMATOKRIT 200 HEMATOLOGY PACKAGE - 24

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PKG. NO.
	24	Capillaries	--	13,000	16,060	200HEMA24

### HAEMATOKRIT 200 HEMATOLOGY BUFFY COAT PACKAGE - 20

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PKG. NO.
	20	Capillaries	--	13,000	16,060	200HEMA20



Learn more at:

[www.hettweb.com/haematokrit200](http://www.hettweb.com/haematokrit200)



Discover more solutions online at: [www.hettweb.com/centrifuges](http://www.hettweb.com/centrifuges)  
**Hettich Lab Technology** | North America | Toll Free: 1 (866) 370 4388 Email: [info@hettweb.com](mailto:info@hettweb.com)

