



EICKEMEYER®
Ultrasonic Doppler

veterinary technology for life
Eickemeyer®

Accurate & reliable

EICKEMEYER® Ultrasonic Doppler

Reliable and accurate non-invasive blood pressure monitor for small animals.

- ▶ Compact hand-held design
- ▶ High sensibility for a clear signal
- ▶ Adjustable sound volume and headphone connection
- ▶ Focused 8 MHz Probe with wide beam for better flow detection in the artery
- ▶ Probe adapted for positioning on the carpal joint or the underside of the tail
- ▶ Rechargeable batteries are not interchangeable
- ▶ Ultrasonic frequency: 8 MHz
- ▶ Operating time: 6 hours
- ▶ Dimensions (in cm): L 17.5 x W 8 x H 2.5
- ▶ Weight: 300 g

217210



217210

Sphygmomanometer

- ▶ Easy one-hand operation
- ▶ Includes two cuffs: 5 cm and 7.5 cm

217309

Blood Pressure Cuff 1.9 cm

- ▶ For small cats

217307

Blood Pressure Cuff 3.3 cm

- ▶ For cats and small dogs

217308

Blood Pressure Cuff 10 cm

- ▶ For large dogs

217311



217308



217309

EICKEMEYER®

Lysbjergvej 6, Hammelev
6500 Vojens
Denmark

T +45 70 20 50 19
info@eickemeyer.dk

TELEFON +45 70 20 50 19

www.eickemeyer.dk