


X-DR Portable Case M WiFi DT

MOBILE SOLUTION FOR EQUINE PRACTITIONERS

W-BLUE BECOMES THE **NEW WE-BLACK**



MORE RELIABLE TECHNOLOGY

Revolutionary new Lossless AED

1-step Auto Sense without loss of x-ray dose

- **With serious X-Ray dose reduction**
- **With much faster image acquisition**
- **With much higher resolution**



Portable handle



MORE DURABLE TECHNOLOGY

Inside AP (Access Point)

Real 2-component design

Shock proof technology

1-meter drop guaranteed by manufacturer

Specialized detector sleeve

- **Supports real portable horse radiology**
- **Unleashed from cables, antenna**
- **Easier detector positioning**

MORE SUITABLE TECHNOLOGY

Award winning

Slim 15 mm Design

1,85 kg Ultralight, 12x 10 "detector with 11 x 9 "active area

Thin 17 mm detector edge for superior equine foot positioning

On-Panel control display

Clip-on batteries

Slim dual battery charger

Allows more suitable x-raying around the horse

- **Improved detector positioning**
- **Added stability**
- **Fits all ISO equine positioning boxes**
- **Clip-on batteries that fit human hand**



MORE PERFORMING TECHNOLOGY

Seriously improved x-ray dose sensitivity!

Enhanced CsI-detector technology

Low noise electronics design

Lossless data transmission

- **Patient dose reduction**
- **Extending your X-Ray range**
- **Fit for lightest portable generator**

MORE ACCURATE TECHNOLOGY

Finally, the equine HR with 76 microns!

Certified for mammographic application.

HR (high resolution) and HT (high transfer) mode in-one

Easy change from HR mode to HT mode

- **More accurate x-ray radiology**
- **Superb diagnostic detail with 76 microns or 6.52 Lp/mm**
- **Keeping fast image transfer and workflow**
- **Easy Switch to HT if you really want to go fast**





Powerful water resistance feature to protect detector

Waterproof Function

EVS 2430W's excellent waterproof design provides protection against blood, sudden rain or other contaminants.

BENEFIT FROM 30 YEARS OF EXPERIENCE AND OUR NETWORK

- Extensive support with project design
- Modular system configuration
- In-house training for technicians
- Installation support if required
- Second level support
- Maintenance and repair

MORE HEALTHY TECHNOLOGY

More detail with less X-Ray dose

More stability with less retakes

- **More time for treatment through less time to X-Ray**
- **Less X-Ray dose to worry about**
- **More range of imaging without increasing X-Ray dose**

LONG BATTERY LIFE THROUGH INTELLIGENT POWER MANAGEMENT



Portable handle

MORE SUSTAINABLE TECHNOLOGY

Power efficient design

Longer battery through lower power consumption

Easy switch from sleep to active mode

- **Better clinics through better technology**
- **Getting more mileage out of batteries before need to recharge**



■ CERTIFICATE: EN ISO 13485:2016

Local Court Stuttgart: HRB 13550
Bank Details, EUR Account: Raiffeisen-Volksbank Ries eG
Bank Details, USD Account: Raiffeisen-Volksbank Ries eG

Tax ID: 90491 / 61321
IBAN: DE93 7206 9329 0000 0113 04
IBAN: DE24 7206 9329 8000 0113 04

VAT-ID: DE147835160
BIC: GENODEF1NOE
BIC: GENODEF1NOE

Chief Officers
Andreas Oppermann, CEO
Bernd Schleicher, CFO