


Introducing iM3[®] Ergo

The new range of *Ergo* dental instruments feature an innovative handle design with soft silicon inserts to provide added grip and comfort during use.

- High quality stainless steel with sharp cutting edges
- Ergonomic designed colour coded handles
- Winged elevators also available extra short for the smallest hands
- Instruments marked for easy identification
- Lightweight yet durable
- Autoclavable to 135 degrees C
- Five-piece set supplied in autoclavable cassette



1. The colour of the ergonomic handle denotes the shape and function of the instrument:




2. The soft silicon grips provide greater comfort and control and are colour coded to help you identify the size of the instrument at a glance. From the Fine Root Tip Luxator up to the larger 6 mm Winged Elevators, the silicon grips all have a consistent colour indicating the tip size as shown below:



Our Extra Short *Ergo* Winged Elevators and Luxators XS give Comfort and Control for even the **Smallest Hands**

Correct fitting hand instruments save time, improve cutting ability, increase tactile sensitivity and overall control resulting in better ergonomics, comfort, technique and safety.

If hand sizes differ greatly within the practice, the recommended solution is to buy two sets in different sizes. This is safe and ergonomically suitable for each user.



Using the convenient hand size guide, you can choose the Ergo Elevator handle that best suits the size of your hand.

Simply place your hand over the illustration below, lining up your index finger and thumb with the **RED** line as shown.

If the width and or length of your hand fits inside the **GREEN** area we recommend the Ergo XS (green handle/purple handle).
If your hand extends into the **BLUE** area, then choose the standard Ergo Elevators (blue handle).



iM3 *Ergo* Winged Elevator XS

Choose *Ergo* XS

Choose Standard *Ergo*

