


Life's better with **HYDRO PHYSIO**[®]
UNDERWATER TREADMILLS



CANINE TREADMILLS AND POOLS


Delivering Canine Care
with Confidence





Hydro Physio Underwater Treadmills and Pools - Delivering Canine Care with Confidence

With more than 20 years of designing and manufacturing canine underwater treadmills and pools, we support Hydrotherapists in delivering enhanced health, recovery and fitness outcomes for their patients.





HYDRO PHYSIO®

Hydro Physio leading the way

Based in Shropshire, England, Hydro Physio have been supplying quality underwater treadmills and pools for over 20 years. In our recently expanded 200,000 Square Foot facility, the design, manufacturing, sales and service teams are all under one roof, facilitating unrivalled customer support.

Innovative design ensures Hydro Physio treadmills have cutting edge features, whilst ISO compliant manufacturing assures a high quality product. The experienced service team offer expert guidance via telephone support, and fully trained service engineers will visit your facility when required.

Over 1,000 Hydro Physio systems have been installed in over 50 countries worldwide.



Thank you to the team at Hydro Physio for their excellent service from ordering to installation. Everyone was friendly and helpful in explaining the process and setting up the treadmill. My spaniel, Looney, very much enjoyed having the guys around to install the treadmill. Thanks again, I am very pleased with the final product!

Ellen Winter
Veterinary Physiotherapist





The Sectors We Serve

Veterinary Community

Provide patients with joint-friendly exercise that improves recovery, pain management and overall quality of life.

Animal Welfare

Help vulnerable animals regain strength and mobility in a safe, controlled environment.



Hydrotherapy

Enhance pain management, mobility and rehabilitation with our underwater treadmills and pools.

MOD and Police

Support working dogs in maintaining exceptional levels of fitness and agility.

Treadmill Specification


Maximum load	90kg (198lbs) without buoyancy
Maximum water depth	0 - 600mm (24")
Speed range	0.2 - 12kph (0.12 - 7.5mph)
Fill / empty times	Approx 2 minutes to max depth
Temperature range	From ambient - 40°C (104°F)
Treadmill dimensions (Excluding display) E	
Excluding ramps	L 2150mm (85") x W 860mm (34") x H 1265mm (50")
Including ramps	L 3435mm (135") x W 860mm (34") x H 1265mm (50")
Complete system weight	
Dry	650kg (1433 lbs)
Full	2150kg (4740lbs)
Water tank capacity	1500 litres (330 gals, 396 US gals)
Max distance between treadmill and plant equipment	8m (314")
Power	240V single phase 50hz 16A (60hz option available)

Pools Specification

MODEL	EXTERNAL MEASURE MM	INTERNAL MEASURE MM	APPROX VOLUME
MAX-1	5400 x 3650 x 1200	4750 x 3000 x 1200	18m3
MAX-2	6650 x 3650 x 1200	6000 x 3000 x 1200	23m3
MAX-3	9150 x 3650 x 1200	8500 x 3000 x 1200	32m3
MAX-4	10400 x 4900 x 1200	9750 x 4250 x 1200	50m3
MAX-5	11650 x 4900 x 1200	1100 x 4250 x 1200	56m3



Life's better with **HYDRO PHYSIO**[®]
UNDERWATER TREADMILLS




..... **Get in Touch**

010 583 50 09

info@eickemeyer.se

eickemeyer.se



OHS 672496



FM 25162

Delivering Canine Care
with Confidence

